


# Patient Slings



**DURABLE**





**Fullback - 1  
Fortrend**

**Page 3**



**Fullback - 2  
Enriched Cotton**

**Page 4**



**Hammock  
High Impact**

**Page 5**



**Toileting (Hygiene)  
Fortrend**

**Page 6**



**Sit 2 Stand  
Easy Applicator**

**Page 7**



**Walking Sling  
Easy Applicator**

**Page 8**



**Walking Pants  
Easy Applicator + High Impact**

**Page 9**



**Lifting, Turning & Repositioning**

**Page 10**



## FULLBACK - 1

Sling Pictured

FORTREND



Sizes  
P S M L XL  
300 Kg  
660 Lbs



Sizes  
XXL  
520 Kg  
1145 Lbs



Breathable  
100%

Also available in:

PATIENT SPECIFIC

HYBRID

EASY APPLICATOR

SUBMARINER

FORTREND

HIGH IMPACT

ENRICHED COTTON

**The Fullback Sling is suitable for most patients who require hoisting.**

The Fullback Sling offers full upper body and head support and stability combined with excellent torso and leg support and stability. The Slings are available in 2 Styles:

- With Straps (Loops)
- With Key Hole Opening Clips

The Fullback Range offers exceptional product features and value including:

- Distribution and dispersion of the patient's weight over the entire sling.
- Minimal stretching under load bearing weight and pressure
- Colour Co-Ordinated Displacement Straps
- Modesty Strap
- Sling Applicator Pocket
- Soft, Turned Webbing Strap Attachment on Sling
- Closed Cell Padded Legs - the padding will not disintegrate or pull apart in laundry
- Suitable for use in Commercial and Institutional Laundries - excellent durability.



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK




[www.rhino-slings.com](http://www.rhino-slings.com)



## FULLBACK - 2

ENRICHED COTTON

- 1 - Permanent Spines - No More Lost Spines!
- 2 - Anti Twist Central Gravity Lift Straps - prevents the sling from twisting and lifts the patient from the centre
- 3 - Sling Applicator Pocket - help with the Insertion and Application of the Sling
- 4 - Pleats to form a Bucket Seat for Safe, Secure Transfers
- 5 - Modesty Straps that also converts into an Amputee Sling
- 6 - Padded Legs using Closed Cell Foam that does not flake or deteriorate in Laundry



Sizes  
P S M L XL  
300 Kg  
1100 Lbs



Sizes  
P S M L XL  
520 Kg  
1100 Lbs



Breathable  
100%

Also available in:

PATIENT SPECIFIC

HYBRID

EASY APPLICATOR

SUBMARINER

FORTREND

HIGH IMPACT

ENRICHED COTTON

Sizes  
S M L XL XXL



ENRICHED COTTON

EASY APPLICATOR



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)



## HAMMOCK

Sling Pictured

HIGH IMPACT



Sizes  
P S M L XL  
300 Kg  
660 Lbs



HYBRID SLING

SUBMARINER SLING

HIGH IMPACT SLING




The Hammock Range offers exceptional product features and value including:

- Distribution and dispersion of the patient's weight over the entire sling.
- Minimal stretching under load bearing weight and pressure
- Colour Co-Ordinated Displacement Straps
- Soft, Turned Webbing Strap Attachment on Sling
- Main Product Label - Full Colour Embroidery, Permanent Colour
- Product Re-Order Label - Permanent Ink - Black & White
- Suitable for use in Commercial and Institutional Laundries - excellent durability.



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)



## TOILETING (Hygiene)

Sling Pictured

FORTREND



Sizes  
P S M L XL  
300 Kg  
660 Lbs



Breakdown  
100%

Also available in:

PATIENT SPECIFIC

HYBRID

EASY APPLICATOR

FORTREND

HIGH IMPACT



**The Toileting (or Hygiene) Sling provides easy access and excellent Patient Stability and Comfort.**

It is used primarily for Toileting (including Cleaning and Hygiene) although can be used as a Transport Sling.

The Sling range is strong and durable and provides the following features:

- Ergonomic Lower Back Design
- Non Slip (Non Latex) Central Back Material
- Fully Padded
- Additional Contoured Tubular Underarm Padding
- A Trinity Sling - Patient Safety and Security System



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)





## SIT 2 STAND

Sling Pictured

EASY APPLICATOR



The Sit 2 Stand Style is designed for use with “Sit-to-Stand” (or “Stand to Pivot”) hoists.

For use with patients who have difficulty in standing up on their own e.g. “stroke” patients, or with single or double arm amputees, and single leg amputees. Can also be used as an aid in physiotherapy to help facilitate standing & walking practice, for example strengthening the legs and feet muscles & promoting patient confidence.

**Please note:** the patients must be able to actively participate in achieving a standing or a semi-standing position. As with all Aids, a risk assessment should be undertaken to determine patient suitability.

The Sling range is strong and durable and provides the following features:

- Ergonomic Lower Back Design
- Non Slip (Non Latex) Central Back Material
- Fully Padded
- Additional Contoured Tubular Underarm Padding
- A Trinity Sling - Patient Safety and Security System



Also available in:

PATIENT SPECIFIC

HYBRID

EASY APPLICATOR

FORTREND

HIGH IMPACT



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)



## WALKING SLING

Sling Pictured

EASY APPLICATOR



The Sit 2 Stand Style is designed for use with “Sit-to-Stand” (or “Stand to Pivot”) hoists.

For use with patients who have difficulty in standing up on their own e.g. “stroke” patients, or with single or double arm amputees, and single leg amputees. Can also be used as an aid in physiotherapy to help facilitate standing & walking practice, for example strengthening the legs and feet muscles & promoting patient confidence.

**Please note:** the patients must be able to actively participate in achieving a standing or a semi-standing position. As with all Aids, a risk assessment should be undertaken to determine patient suitability.

The Sling range is strong and durable and provides the following features:

- Ergonomic Lower Back Design
- Non Slip (Non Latex) Central Back Material
- Fully Padded
- Additional Contoured Tubular Underarm Padding
- A Trinity Sling - Patient Safety and Security System



Also available in:

HYBRID


EASY APPLICATOR

FORTREND

HIGH IMPACT



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)



## WALKING PANTS

Sling Pictured

EASY APPLICATOR  
+  
HIGH IMPACT



Also available in:

HYBRID  
EASY APPLICATOR  
FORTREND  
HIGH IMPACT

Sizes  
S M L XL

**The Walking Pants are designed as an aid to assist with patients who have difficulty in walking.**

The Walking Pants may also be used in rehabilitation and gait training.

The Walking Pants are generally designed for use with Ceiling hoists with "Tracks" because of Height Restrictions - "Floor" (Mobile) Lifter/Hoist may not be suitable for Tall Patients. Ceiling hoists with Tracks also offer more versatility and freedom of movement.

Please note that the Pants are intended for use with patients who are able to actively participate in achieving a standing position and are able to support their weight.

The Sling range is strong and durable and provides the following features:

- Anatomically Shaped for better patient comfort
- Fully Padded
- Easy to apply




Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)



Rhino Consultant and Facilitators Ltd  
10 Gregston Trade Centre  
Birmingham Road  
Oldbury - B69 4EX UK



[www.rhino-slings.com](http://www.rhino-slings.com)